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4. CARY A.BROWN. Occupational Therapy International, 9(2), 121-130, 2002.

# SCARLESS

GEL

**A silicone gel  
to help flatten and soften scars**

**Indication:** For the prevention and treatment of scars caused by any trauma to the skin (e.g.: injuries, burns, wounds, plastic and cosmetic surgery, animal bites, etc.) SCARLESS flattens the edges of the scar, softens it and improves its color. SCARLESS is also effective in flattening the edges and improving the color of scars caused by acne (pimples of puberty).

**Composition:** Bio-inert silicone ingredients, such as: silicone dioxide and polysiloxane

**What is a scar?** A scar is a mark formed as a result of the biological process of wound healing in the skin and other tissues of the body. A skin scar is caused when the deepest layer of the skin, the dermis, is damaged. In order to repair the damage caused, the body produces new collagen fibers. This process leads to a scar. Since the body cannot reconstruct the tissue as it was, the new tissue will have a texture and quality different from the surrounding tissue.

There are two types of scars:

A hypertrophic scar - A raised red lump on the skin, which does not grow beyond the local wound.

A keloid scar - A serious condition where the body produces collagen indefinitely, resulting in a large, tumorous (although benign) lump.

A scar may also take the form of a sunken recess in the skin. Such recesses are formed when underlying structures supporting the skin, such as fat or muscle tissues, are lost<sup>(1)</sup>.

**About the silicone gel in SCARLESS:** SCARLESS contains a mixture of medical silicones designed for application on the scars. The silicone gel forms a selective barrier separating the skin surface from its surroundings, thus creating a breathable barrier, impermeable to water, that blocks the passage of bacteria. This enables the immune system to repair the damage more effectively and allows normal activity and collagen production<sup>(2)</sup>. In addition, the silicone gel moisturizes the scar, and as a result the scar heals more quickly, and has a softer, more delicate and esthetic appearance.

**Dosage and method of use:** Begin using SCARLESS only after the wound appears closed, or after the sutures have been removed. Do not apply SCARLESS on mucous membranes or on other preparations.

Safe for use in adults and children. For external use only. Clean and dry the scar. Apply a thin layer of the gel on clean and dry skin every morning and evening. Wait about 10-15 minutes for the gel to dry. After the gel has dried, you may get dressed or apply make-up to the treated area.

**Additional guidelines:** It is recommended to begin using SCARLESS shortly after wound closure in order to achieve better results. It may also be used for old scars<sup>(3)</sup>; for scars four or more years old, the results will be negligible.

The results of using SCARLESS will be apparent about one month after beginning treatment. For optimal results, the gel should be used for at least three months.

**The benefits of SCARLESS gel:**

- Easy and convenient to use; the gel dries rapidly after application.
- Transparent, colorless and odorless, thus does not interfere with your daily life.
- Also suitable for use on children<sup>(4)</sup> and in people with sensitive skin.

**Storage conditions:** Store at room temperature (25°C). Keep out of the reach of children. Read the patient leaflet before use.

**Medical Device registration number:** 14380403

**Manufactured by:** MediCure Technology Ltd.,  
5 Ha-Shita St., Industrial Zone, Caesarea.

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silicone dioxide and polysiloxane

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